



检测报告

Test Report

送检企业: 四川金澜科技有限公司

样品名称: 超纯水设备产水水样

检验类别: 委托检验 (在线检测)

厦门泓益检测有限公司

Xiamen Hongyi Testing Technology Co.Ltd

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of Hongyi Testing, this test report shall not be copied except in full and published as advertisement Hongyi Physical & Chemical Lab.

厦门泓益检测有限公司
Xiamen Hongyi Testing Co., Ltd.

地址: 福建省厦门市同安区轻工食品工业区美禾路 99 号

ADDRESS: No. 99 Meihe Road, Light Industry Food Industrial Zone, Tongan District, Xiamen City, Fujian Province,

邮箱: hongyotest@123.com



报告编号: HYI20260821A1609

样品名称	:	超纯水设备产水水样	规格型号	:	/
商标	:	金蓝超纯	样品数量	:	2 PCS
委托单位	:	四川金澜科技有限公司			
委托单位地址	:	成都市成华区龙潭工业园成致路 48 号 3 号楼 2 楼			
生产单位	:	四川金澜科技有限公司			
生产单位地址	:	成都市成华区龙潭工业园成致路 48 号 3 号楼 2 楼			
检验类别	:	委托检验（在线检测）	测试日期	:	2026 年 08 月 10 日-2026 年 08 月 19 日
签发日期	:	2026 年 08 月 21 日			
检验依据	:	1、GB/T 6682-2008《分析实验室用水规格和试验方法》 2、GB/T 33087-2016《仪器分析用高纯水规格及试验方法》			
检验项目	:	见后续			
样品特征和状态	:	完好			
环境温度	:	24.6℃	环境湿度	:	63.8%
测试结果	:	本次委托检验，所检项目全部符合标准要求			
样品描述	:	本次检测水样来自于 JL-UPH-IV-20T 设备产出水样（主测）；试验前超纯水、纯水水样样品完好。 附测：JL-UPH-IV-5T、JL-UPH-IV-10T、JL-UPH-IV-10TM、JL-UPH-IV-15T、JL-UPH-IV-20TM、 JL-UPH-IV-30T、JL-UPH-IV-40T、JL-UPH-IV-40L、JL-UPH-IV-50L、JL-UPH-IV-50LM、JL-UPH-IV-60L、 JL-UPH-IV-70L、JL-UPH-IV-80L、JL-UPH-IV-100L、JL-UPH-IV-100LM、JL-UPH-IV-120L、 JL-UPH-IV-120LM、JL-UPH-IV-150L、JL-UPH-IV-180L、JL-UPH-IV-200L、JL-UPH-IV-200LM、 JL-UPH-IV-250L、JL-UPH-IV-250LM、JL-UPH-IV-300L、JL-UPH-IV-300LM			

报告编制: 张华伟

报告审核: 杨华群

报告签发: 杨华群

检测专用章

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of Hongyi Testing, this test report shall not be copied except in full and published as advertisement Hongyi Physical & Chemical Lab.



报告编号: HYI20260821A1609

序号	检测项目	检测指标	检测结果	单项判定
一	一级水 (超纯水)			合格
1	电阻率 (25℃)	$\geq 18 \text{ M}\Omega \cdot \text{cm}$	18.25	合格
2	电导率 (25℃)	$\leq 0.1 \mu\text{S/cm}$	0.055	合格
3	吸光度 (254nm, 1cm 光程)	≤ 0.001	0.0005	合格
4	总有机碳 (TOC)	$\leq 50\text{ppb}$	2.62	合格
5	钠离子 (Na^+)	$\leq 1\text{ppb}$	0.01	合格
6	氯离子 (Cl^-)	$\leq 1\text{ppb}$	0.01	合格
7	硅 (Si)	$\leq 10\text{ppb}$	1.4	合格
8	微量元素 (Mn/Fe/Cu/Zn)	$< 0.1\text{ppb}$	0.02	合格
9	细菌总数	$< 1 \text{ CFU/mL}$	未检出	合格
10	颗粒物 ($>0.2 \mu\text{m}$):	$< 1 \text{ 个/mL}$	< 1	合格
11	微粒及细菌去除率	$> 99\%$	99.8%	合格
12	内毒素	$< 0.03 \text{ EU/mL}$	0.01	合格
13	核糖核酸酶	$< 1 \text{ pg/mL}$	未检出	合格
14	脱氧核糖核酸酶	$< 5 \text{ pg/mL}$	3	合格
二	三级水 (纯水)			合格
1	电阻率 (25℃)	$\geq 0.2 \text{ M}\Omega \cdot \text{cm}$	0.46	合格
2	电导率 (25℃)	$\leq 5 \mu\text{S/cm}$	2.16	合格
3	可氧化物质 (O 计, mg/L)	≤ 0.4	0.3	合格
4	105℃蒸发残渣 (mg/L)	≤ 2.0	1.0	合格
5	pH 值 (25℃)	5.0~7.5	6.5	合格

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of Hongyi Testing, this test report shall not be copied except in full and published as advertisement Hongyi Physical & Chemical Lab.

厦门泓益检测有限公司

Xiamen Hongyi Testing Co., Ltd.

地址: 福建省厦门市同安区轻工食品工业美禾路 99 号

ADDRESS: No. 99 Meihe Road, Light Industry Food Industrial Zone, Tongan District, Xiamen City, Fujian Province,

邮箱: hongyotest@123.com



样品照片



JL-UPH-IV-20T

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of Hongyi Testing, this test report shall not be copied except in full and published as advertisement Hongyi Physical & Chemical Lab.

厦门泓益检测有限公司
Xiamen Hongyi Testing Co., Ltd.
地址: 福建省厦门市同安区轻工食品工业区美禾路 99 号
ADDRESS: No. 99 Meihe Road, Light Industry Food Industrial Zone, Tongan District, Xiamen City, Fujian Province,
邮箱: hongyotest@123.com



报告编号: HYI20260821A1609

注 意 事 项:

1. 报告无“检验报告章”或检验单位公章无效。
2. 未经本实验室书面同意,不得部分复制本报告。
3. 报告无主检、审核、批准人签字无效。
4. 报告涂改无效。
5. 对检测报告若有异议,应于收到报告之日起十五日内向检测单位提出,逾期不予受理。
6. 委托检测仪对来样负责。

地址:福建省厦门市同安区轻工食品工业区美禾路 99 号

Add: No. 99 Meihe Road, Light Industry Food Industrial Zone, Tongan District, Xiamen City, Fujian Province,

报告结束

This report shall not be altered, increased or deleted. The results shown in this test report refer only to the sample(s) tested. Without written approval of Hongyi Testing, this test report shall not be copied except in full and published as advertisement Hongyi Physical & Chemical Lab.

厦门泓益检测有限公司
Xiamen Hongyi Testing Co., Ltd.

地址:福建省厦门市同安区轻工食品工业区美禾路 99 号

ADDRESS: No. 99 Meihe Road, Light Industry Food Industrial Zone, Tongan District, Xiamen City, Fujian Province,

邮箱: hongyotest@123.com